

I AM
CONFIDENT




SOPIX SERIES

**A successful X-ray
every time with minimal
exposure to radiation**



STRIKING CONTRAST FOR A MORE RELIABLE DIAGNOSIS





MORE INVENTIVE Better differentiation of the dental tissue

SOPIX® sensors surpass the limits of radiological examinations by offering **greater differentiation of dental tissue.**

This technological achievement is called **FIBER2PIXEL®.**

WITH FIBRE - FIBER2PIXEL®

WITHOUT FIBRE



2 FIBER PIXEL

Differentiation of the dental tissue

FIBER2PIXEL® technology is based on the use of **broad spectrum optical microfibres** for the guided transmission of photon emissions in order to provide **highly contrasted images.**

LESS INVASIVE

A more reliable diagnosis

The different tooth anatomic structures, such as the bone, roots, pulp... are highlighted with **extreme precision** on the image.

Your diagnosis is **faster** and **more accurate!**



THE PERFECT FIT TO YOUR CLINICAL APPLICATIONS

Endodontics



Pedodontics



Cariology




Periodontics




Periapical



Bitewing



Implantology



HIGH-QUALITY IMAGES

With **FIBER2PIXEL**® technology, SOPIX® sensors provide **accurate images** and **striking contrast** to ensure a **reliable diagnosis**.

DESIGNED FOR YOUR PRACTICE

Two sizes are available depending on **patient morphology** and **clinical applications**.

Scale 1



Size 2



Size 1



ACTEON IMAGING SUITE

THE NEW GENERATION OF IMAGING SOFTWARE

Extremely user-friendly, ACTEON Imaging Suite software offers **advanced X-ray image processing tools**.

ACTEON Imaging Suite is delivered with each SOPIX and is compatible Windows® and Mac®.


A QUALITY IMAGE EVERYTIME WITH MINIMAL EXPOSURE TO RADIATION

CUTTING EDGE TECHNOLOGY

Available in all SOPIX® series sensors, patented **Ace** technology (Automatic control exposure) analyses in real-time, the amount of X-rays accumulated by the sensor. It automatically freezes the image acquisition as soon as the sensor receives the radiation required to produce the perfect image.

Eliminate the risk of over exposing the image!

Combined with the X-Mind® unity intraoral X-ray generator, SOPIX inside with ACE technology **limits the emission of x-rays** during the acquisition to the necessary amount for the patient's morphology. It uses the **minimum dose** required to provide a high-quality image.



"ACE is the combination of advanced sensor technology, digital power electronics and the know-how of two diagnostic imaging divisions. The synergy between La Ciotat (FRANCE) and Milan (ITALY) R&D teams gave birth to an innovative concept focused on patients, with outstanding image quality."



FOR A SAFER PROCESS

With SOPIX Series sensors and its patented ACE technology, you acquire **successful X-rays every time**, meaning reliable and accurate diagnosis. You **save time** avoiding the need for retakes.

Whilst using X-Mind unity intraoral X-ray generator with SOPIX inside, the patients **receive the minimum required dose for their dental morphology**. You protect your patients and your staff from unnecessary radiation.

PATIENT AND STAFF

OPTIMAL PROTECTION



STOP EXCESSIVE RADIATION

The communication between the X-Mind unity and SOPIX inside sensor provides **unique benefits**.


When SOPIX inside has received enough energy to provide an **exceptional image quality**, it tells the X-Mind unity to **stop the X-ray emission**.



Effective protection for minimal exposure

The patient only receives the necessary dose adapted for their dental morphology, which **protects them from unnecessary exposure**.

ACTEON Imaging Suite, always one step ahead



ACTEON Imaging Suite systematically records the **X-Mind unity settings** as well as the **effective dose received by the patient** for each acquisition.

This ensures **permanent traceability** for every patient.



EXCLUSIVE TRACEABILITY

Outstanding working comfort

Through direct integration of SOPIX inside into X-Mind unity, **connecting cables are hidden** inside the X-ray unit.

The holder places the sensor **safely at easy reach** to prevent it from falling onto the floor.

Your working environment is therefore **more ergonomic and productive**.



OUTSTANDING PERFORMANCE

SMART DESIGN FOR BETTER COMFORT

White side stripes ensure high visibility of the sensor in the dark area of the mouth, to correctly position the X-ray tube perpendicular to the sensor.



Rounded edges and corners for improved patient comfort.

FAST AND EASY TO USE

Save time with a sensor that is **always ready to acquire**.
The image is **displayed immediately**.

NO MORE OVEREXPOSED IMAGES

Available on all SOPIX series sensors,
ACE technology freezes the image during acquisition to protect it from over-exposure.

Acquire perfect image the first time and every time!



THE SOPIX SERIES

SOPIX

With proven quality and reliability, SOPIX produces a high quality image at an affordable price.

The most economic solution of the SOPIX series



SOPIX²

This sensor provides an exceptional image quality, using the most advanced technology.

The solution for optimal performance



This sensor is directly integrated into the X-Mind unity intraoral X-ray generator, resulting in a reduction of X-ray emissions.

The patient's well being is the highest priority



TECHNICAL SPECIFICATIONS

Size 1

External dimensions 25 x 39mm
 Active surface area 600mm² (20 x 30mm)
 Number of pixels 1.50million

SOPIX / SOPIX inside system

Technology CMOS + scintillator+ optic fibre
 Pixel size 20µm x 20µm
 Theoretical resolution 25lp/mm
 Real resolution >12lp/mm
 Supplied imaging software ACTEON Imaging Suite
 TWAIN module Yes

SOPIX / SOPIX² USB connection

Connection USB 2.0
 Total cable length 3.70m

Windows® minimum configuration required

Operating system Windows® 7
 Processor Quadcore 2.6 Ghz
 RAM 4 GB
 Hard disk 300 GB
 USB ports 2 USB 2.0 Hi-Speed ports
 Graphic card OpenGL 2.1 or better
 alternatively DirectX 9 or 11 Graphics Device
 USB Chipset Intel® or NEC® / RENESAS®
 Screen resolution 1600 x 1024
 Ethernet board 100 Mbps - 1 Gbps

Mac® minimum configuration required

Computer MacBook® Pro 13.3" or iMac® 21.5"
 Operating system 10.12 Sierra
 Processor Quadcore 2.6 Ghz
 RAM 4 GB
 Ethernet board 1 Gbps

Size 2

External dimensions 31 x 42mm
 Active surface area 884mm² (26 x 34mm)
 Number of pixels 2.21 millions

SOPIX² / SOPIX² inside system

Technology CMOS + scintillator + optic fibre
 Pixel size 20µm x 20µm
 Theoretical resolution 25lp/mm
 Real resolution >18lp/mm
 Supplied imaging software ACTEON Imaging Suite
 TWAIN module Yes

SOPIX inside / SOPIX² inside USB connection

Connection USB 2.0
 Sensor cable length 0.70m

Windows® recommended configuration

Operating system Windows® 10
 Processor Quadcore 2.6 Ghz
 Ram 8 GB
 Hard disk 1 TB
 USB ports 4 USB2 Hi-Speed ports
 Graphic card Dedicated graphics card with
 a t least 1 GiB memory
 USB Chipset Intel® or NEC® / RENESAS®
 Screen resolution 1920 x 1080 for optimal planning
 or better
 Ethernet board 1 Gbps

Mac® recommended configuration

Computer iMac® 27"
 Operating system 10.14 Mojave
 Processor Quadcore 2.6 Ghz
 RAM 8 GB
 Ethernet board 1 Gbps

For Yosemite and El Capitan operating systems, a Mac computer from 2013 or later is required.

Note: In the case of SOPIX inside and SOPIX² inside, the IEC 60601-2-65 norm requires for each X-Ray intraoral system with an onboard digital sensor to use a square collimator.

The medical devices for dental care SOPIX Series are of class IIa and manufactured by SOPRO, notified body GMED CE0459, X-Mind unity and AIS are class IIb and manufactured by DE GÖTZEN, notified body IMQ - CE 0051. These medical devices are not refunded by health insurance organizations. Read carefully the instructions on the labelling before use.

SOPIX®, X-Mind®, FIBER2PIXEL® and SOPRO® are registered trademarks of SOPRO.
 X-Mind® is registered trademarks of DE GÖTZEN.
 "All other trademarks cited herein are the property of their respective owners"



www.acteongroup.com

For more information, please contact:

SOPRO S.A. | A company of ACTEON Group
ZAC Athélia IV | Avenue des Genévriers | 13705 LA CIOTAT cedex | FRANCE
Tél + 33 (0) 442 98 01 01 | Fax + 33 (0) 442 71 76 90
E-mail: info@sopro.acteongroup.com | www.acteongroup.com

